

ZINC PREP.

21. 11. 85

11

METHOD SHEET FOR PASTEURISED FULL PILOT-SCALE PRODUCTION TESTBATCH ZHTDAY 1

1. Tris extract (pH 7.0) 830 ml l *
2. Titrate slowly from 7.0 to 6.7 6.7 with 0.05N HCl
Use E₁ vibromixer
Volume used — ML's
Total volume 830 ML's (V₁)
3. Zinc/heparin precipitation, pH 6.7
Use E₁ vibromixer
~~0.04 ML, Heparin/l(V₁) — ML's~~
333.3 ML's, 4mM ZnAc/l(V₁) 277 ML's add slowly.
Final pH, 6.7
4. Equilibrate 10 minutes R.T. with gentle mixing.
5. Centrifuge 10 minutes R.T. 4000 rpm.
Post-centrifugation, zinc-supernatant 980 ml (V₂)*
6. Mixture of citrate and calcium chloride to 0.02M and 2.5 mM
a) $41.7 \text{ MLS/l(V}_2) = 41.7 \times \frac{980}{1000} \text{ l(V}_2) = 40.9 \text{ MLS 0.5M Citrate}$
b) $11.1 \text{ MLS/l(V}_2) = \frac{11.1}{5.05} \times \frac{980}{1000} \text{ l(V}_2) = 4.9 \text{ MLS 0.5M Calcium Chloride.}$
Mix 40.9 MLS 0.5M CaCl₂ with 4.9 MLS 0.5M Citrate
Add slowly to (V₂) use E₂ vibromixer
Final volume 1030 l *
- Adjust to pH 6.7 and 0.45 filter 2x20ml aliquots for rest in 1 lump.
7. ~~Clarify through a-NBP cartridge — l*(V₃)~~ 800 ml
8. ~~Addition of glycine and sorbitol stabilisers~~
 - a) $50 \text{ g/l(V}_3) = 50 \text{ g} \times \text{— l(V}_3) = \text{— g glycine.}$
Add slowly mix with E₂ vibromixer.
Keep pH above 7.0 during glycine addition with 1N NaOH and adjust to final pH 7.5.
 - b) $1850 \text{ g/l(V}_3) = 1850 \text{ g} \times \text{— l(V}_3) = \text{— g sorbitol}$
Add over 30 minutes with stirrer at 70 rpm. INCREASE temperature GRADUALLY from +20°C to +40°C.
Fill glycine/sorbitol VIII (G+S) into MSE bottles
— x — ML'S

SNBTS DOCUMENT REQUEST No:

2010/00045

Zinc Prep. 21.11.85.

	VOL ML'S.	1-ST FVIII U/ML	u	2-ST FVIII U/ML	u.	TPg/2g	TPg/2g	g	% FVIII RECOVERY	1-ST. % FVIII RECOVER
TS ₁ PH7.0	830	8.96	7437	7.04	5843	34.8 (29)	20.5 (17)			
						(28884)				
TS ₁ PH7.467	830	9.14	7586	7.23	6001	35.3 (29)	17.9			
Z-S	980	5.94	5821	4.97	4871	16.5 (16)	4.9		77	81
Z-S +cit/ser.	1030	5.79	5964	5.26	5418	17.3 (17.8)	7.2		79	90
Post 1.2μ	-	5.70		5.43		16.7	7.0			
" 0.8μ	-	5.34		4.64		16.6	7.0			
" 0.45μ	860	5.33	4584	4.82	4145	16.8 (14.5)	7.1 (6.1)		60	70

% loss TP.
50% TP
64

METHOD SHEET FOR PASTEURISED EVIII PILOT-SCALE PRODUCTION TEST

BATCH ZHT

DY 5069

2.12.85

DAY 1

1. Tris extract (pH 7.0) 525 ml l *
2. Titrate slowly from 7.1 to 6.7 6.7 with 0.05N HCl
Use E, vibromixer
Volume used _____ ML's
Total volume 500 ML's (V₁)
3. Zinc/heparin precipitation, pH 6.7
Use E, vibromixer
0.04ML, Heparin/l(V₁) - ML's
333.3 ML's, 4mM ZnAc/l(V₁) 167 ML's add slowly.
Final pH, 6.65
4. Equilibrate 10 minutes R.T. with gentle mixing.
5. Centrifuge 10 minutes R.T. 4000 rpm.
Post-centrifugation, zinc-supernatant 640 l(V₂)* (620)
6. Mixture of citrate and calcium chloride to 0.02M and ^{2.5}~~5.5~~ mM
a) $41.7 \text{ MLS/l(V}_2) = 41.7 \times \text{l(V}_2) = \underline{25.9} \text{ MLS } 0.5\text{M Citrate}$
b) $\frac{11.1}{5.5} \text{ MLS/l(V}_2) = \frac{11.1}{5.5} \times \text{(V}_2) = \underline{3.4} \text{ MLS } 0.5\text{M Calcium Chloride.}$
Mix 3.4 MLS 0.5M CaCl₂ with 26.0 MLS 0.5M Citrate
Add slowly to (V₂) use E2 vibromixer adjust pH to 7.4
Final volume 670 l *
7. Clarify through a NBP cartridge filtered no problem ^{0.8} 0.45 μ post-filtr vol 630ml l*(V₃)
8. Addition of glycine and sorbitol stabilisers
a) $50 \text{ g/l(V}_3) = 50 \text{ g} \times \text{l(V}_3) = \text{_____ g glycine.}$
Add slowly mix with E2 vibromixer.
Keep pH above 7.0 during glycine addition with 1N NaOH and adjust to final pH 7.5.
b) $1850 \text{ g/l(V}_3) = 1850 \text{ g} \times \text{l(V}_3) = \text{_____ g sorbitol}$
Add over 30 minutes with stirrer at 70 rpm. INCREASE temperature GRADUALLY from +20°C to +40°C.
Fill glycine/sorbitol VIII (G+S) into MSE bottles
_____ x _____ ML'S

NYS069.

2.12.85.

Vol	1ST FVIII μ ML	μ	2ST FVIII μ ML	μ	TP g/L	g	TP g/L	g	%RECOVERY 2-ST VIII	%RECOVERY TP 1-ST VIII
TB ₁ pH7.0 (S2S)	2.8	4095	7.1	3728	29.7	106	14.7	77		
TB ₁ pH6.7 (S2S)	7.1	3722	7.2	3780	24.1	127	11.2	5.9		
Z-S (640)	4.4	2216	4.5	2280	11.3	72	2.5	16		
Z-S+ClT+Ca (670)	4.0	2680	4.4	2942	11.2	75	3.1	2.1		
pH7.4 " (670)	4.9	2223	4.7	3149	11.1	74	2.9	1.9	88.0	83.0
POST-0.8 μ FILT.	4.0		4.1		10.1		2.5			
POST-0.45 μ (630) FILT.	3.72	2344	3.7	2331	9.8	62	2.8	1.2	63.0	62.0

	%RECOVERY 2-ST VIII	%RECOVERY 1-ST VIII	%LOSS TP	%LOSS Tq.
POST Z-S+ClT+Ca	88	83	41	68
POST-0.45 μ FILT.	63	62	51	69