

ENDOTOXIN TREATED FACTOR VIII IN VIVO RECOVERY AND * LIFE STUDIES HAEMOPHILIA

PATIENT'S NAME..

HOSPITAL.....R.I.E.

D.O.B.....

BODY WEIGHT.....68

Kg. HEIGHT.....5' 9"

CONSULTANT.....

DR. L. L. L. L. L.

PARTICULARS OF INFUSION:

Batch No.....024100631

No. of Vials.....10 Volume infused.....227 mls. Total dose.....2,500 u

If local VIII C assay on infusate was performed please give value and also value & source of standard.....

Solubility time.....35 minutes. Infusion time.....36 minutes

	1 HOUR 2 HOUR 3 HOUR 4 HOUR 5 HOUR 6 HOUR 7 HOUR 8 HOUR 9 HOUR 10 HOURS									
	Before	10 min.	20 min.	30 min.	60 min.	120 min.	180 min.	240 min.	300 min.	360 min.
TEMPERATURE	36.7	36.7	36.7	36.2	36.7	36.7	36.7	36.0	36.5	36.4
BLOOD PRESSURE	145/80	145/85	130/80	128/80	130/80	130/80	130/80	125/86	125/86	135/90
PULSE	87	80	76	80	79	84	92	95	88	91
Factor VIII C	0.11				0.83			0.61	0.61	0.25
Factor VIII Cag										
Factor VIII Rag										
Anti H3 (litre)										
10 mls plasma stored at -30°C										

NOTES: (1) Dose should be infused in 20 minutes.

(2) 10 mls plasma should be aliquoted (1 ml) and stored for future studies.

(3) Patients selected for this study should be haemostatically stable.

(4) PLEASE RETURN COMPLETED FORM TO DR F BOULTON, EDINBURGH & SOUTH EAST BLOOD TRANSFUSION SERVICE, ROYAL INFIRMARY, EDINBURGH.

Clinical Comments

pin point results on 10 mls plasma stored at -30°C