

NY 776

pkch 11/11/85.

Unheated

- 1 NY 776
- 2 " + 200 mM lysine
- 3 " + 2% sucrose

Heated (80°C / 72 hr)

- 3 NY 776
- 4 " + 2% sucrose
- 5 " + 5% sucrose
- 6 NY 776 + 100 mM lysine
- 7 " 200
- 8 " 300
- 9 NY 776 + 100 mM lysine + 2% sucrose
- 10 " 200 " "
- 11 " 300 " "
- 12 NY 776 + 5% sucrose + 100 mM lysine
- 13 " 200
- 14 " 300
- 16 NY 776 + 10% sucrose + 100 mM lysine.
- 17 " 200
- 18 " 300
- 19 } NY 195, lysine, SPS
- 20 }